

Protocol

Chromatin Immunoprecipitation

James DeCaprio and Thomas O. Kohl

Chromatin immunoprecipitation, commonly referred to as ChIP, is a powerful technique for the evaluation of *in vivo* interactions of proteins with specific regions of genomic DNA. Formaldehyde is used in this technique to cross-link proteins to DNA *in vivo*, followed by the extraction of chromatin from cross-linked cells and tissues. Harvested chromatin is sheared and subsequently used in an immunoprecipitation incorporating antibodies specific to protein(s) of interest and thus coprecipitating and enriching the cross-linked, protein-associated DNA. The cross-linking process can be reversed, and protein-bound DNA fragments of optimal length ranging from 200 to 1000 base pairs (bp) can subsequently be purified and measured or sequenced by numerous analytical methods. In this protocol, two different fixation methods are described in detail. The first involves the standard fixation of cells and tissue by formaldehyde if the target antigen is highly abundant. The dual cross-linking procedure presented at the end includes an additional preformaldehyde cross-linking step and can be especially useful when the target protein is in low abundance or if it is indirectly associated with chromatin DNA through another protein.

MATERIALS

It is essential that you consult the appropriate Material Safety Data Sheets and your institution's Environmental Health and Safety Office for proper handling of equipment and hazardous materials used in this protocol.

RECIPES: Please see the end of this protocol for recipes indicated by <R>. Additional recipes can be found online at <http://cshprotocols.cshlp.org/site/recipes>.

Reagents

Agarose gel (1%)
Antibody (target dependent)
Block solution (PBS containing 0.5% BSA)
Cells, 10^8 healthy viable cells
 Monolayer-cultured cells: $10 \times 15 \text{ cm}^2$ plates
 Suspension-cultured cells: 0.5×10^6 to 1.0×10^6 cells/mL
Cells, monolayer cultured ($5 \times 15\text{-cm}^2$ plates) (for dual cross-linking procedure only)
Disuccinimidyl glutarate (DSG; 0.5 M) <R> (for dual cross-linking procedure only)
DNA loading buffer (6×) <R>
EBC lysis buffer for ChIP <R>
Elution buffer for ChIP <R>

From the Antibodies collection, edited by Edward A. Greenfield.

© 2020 Cold Spring Harbor Laboratory Press

Cite this protocol as *Cold Spring Harb Protoc*; doi:10.1101/pdb.prot098665



Formaldehyde (37%)
Formaldehyde solution (11%) <R>
Glycine (2.5 M)
N-Lauroyl-sarcosine (10% stock solution)
Liquid nitrogen (LN₂)
Lysis Buffer 1 for ChIP <R>
Lysis Buffer 2 for ChIP <R>
Lysis Buffer 3 for ChIP <R>
PBS (pH 7.2)
PMSF (for dual cross-linking procedure only)
Protease inhibitors (Set I, Calbiochem catalog no. 539131; or similar)
Proteinase K (Roche)
RIPA buffer solution for ChIP <R>
RNase A (Sigma-Aldrich)
Salmon sperm DNA (sonicated)
T₅₀E₁₀ (without SDS)
T₅₀E₁₀S₁ mixture
TE buffer containing 50 mM NaCl

Equipment

Aspirator
Cell scraper
Conical tube (50 mL; Falcon)
Dynabeads and magnetic stand
Eppendorf tubes (1.5 mL)
Gel electrophoresis equipment
Humidified incubator set at 37°C, 5% CO₂
Hybridization oven
Magnetic racks
Microcentrifuge
NanoDrop Spectrophotometer
PCR apparatus
PCR purification columns (QIAGEN catalog no. 28104)
Rotating platform
Rotating wheel (for dual cross-linking procedure only)
Sonicator
Vortex

METHOD

Chromatin Immunoprecipitation

This protocol details individual steps sequentially performed during the chromatin immunoprecipitation procedure and is ideally suited for the formaldehyde fixation of chromatin to highly abundant proteins. The initial steps of chromatin preparation and evaluation before the performance of the immunoprecipitation of the target antigen(s) are described in this protocol, followed by step-by-step guidance in the washing and elution of target antigen(s) and the reverse cross-linking of the protein-associated DNA fragments for purification, amplification, and sequencing purposes.



Cross-Linking of DNA and Protein

For Monolayer-Cultured Cells

1. Remove the growth medium from ten 15-cm² cell culture dishes of proliferating cells totaling $\sim 10^8$ cells, and wash the adherent cells once with PBS before aspiration.
Avoid the use of overconfluent or dead cells.
2. Fix the cells by addition of 10 mL of 11% Formaldehyde solution, and gently incubate the cell culture dishes for 10 min at room temperature using a rotating platform.
3. Stop the cross-linking reaction by adding 0.5 mL of 2.5 M glycine to a final concentration of 125 mM. Incubate the cell culture dishes for 5 min at room temperature on a rotating platform, and aspirate the solution completely from each dish.
4. Add 5 mL of ice-cold PBS to individual dishes, and scrape the cells from all plates into two 50-mL conical tubes. Collect cells by centrifugation at 400g for 5 min at 4°C, and aspirate the supernatant.
5. Wash the fixed cell pellet with 5 mL of ice-cold PBS, and pellet cells by centrifugation at 400g for 5 min at 4°C. Aspirate the PBS completely.
6. Proceed to Step 14 or snap-freeze the cells in LN₂ for storage at –80°C.

For Suspension-Cultured Cells

7. Collect cells from a 100-mL culture (0.5×10^6 to 1.0×10^6 cells/mL) by centrifugation at 400g for 5 min, and aspirate the growth medium.
8. Wash the cells once with PBS and pellet by centrifugation at 400g for 5 min. Aspirate the Wash buffer.
9. Resuspend the cell pellet in 20 mL of 11% Formaldehyde solution and gently incubate for 15 min at room temperature using a rotating platform.
10. Stop the cross-linking reaction by adding 1.0 mL of 2.5 M glycine to a final concentration of 125 mM, and incubate the cells for 5 min at room temperature on the rotating platform.
11. Collect the cells by centrifugation at 400g for 5 min at 4°C, and aspirate the supernatant.
12. Wash the fixed cell pellet with 5 mL of ice-cold PBS, and pellet the cells by centrifugation at 400g for 5 min at 4°C before complete aspiration of the PBS.
13. Proceed to Step 14 or snap-freeze the cells in LN₂ for storage at –80°C.

Chromatin Preparation

14. Resuspend cell pellets in 10 mL of Lysis Buffer 1 containing protease inhibitors and incubate for 10 min at 4°C using a rotating platform.
15. Collect cells by centrifuging at 400g for 5 min at 4°C, and resuspend the pellet in 10 mL of Lysis Buffer 2 containing protease inhibitors. Incubate the tube for 10 min at room temperature on a rotating platform.
16. Pellet the cells by centrifuging at 400g for 5 min at 4°C, and resuspend the pellet in 3–4 mL of Lysis Buffer 3 containing protease inhibitors.
17. Sonicate the cells on ice using a 20% duty cycle for 2 min followed by a 1-min cooling period. Repeat this 10–15 times.
The solution will start to appear clear when looking through the tube. The goal is to shear chromatin DNA into fragments ranging 100–600 bp in length. Although chromatin from some cell lines is easy to shear into this range, the shearing of chromatin from other lines is extremely difficult and requires more extensive sonication at higher power settings over more cycles.
18. Add 150 µL of *N*-lauroyl-sarcosine (10% stock solution) and incubate for 10 min at room temperature on a rotating platform.

19. Distribute the mixture into three 1.5-mL Eppendorf tubes and centrifuge at 15,000g for 5 min at 4°C using a microcentrifuge. Collect the supernatant containing the chromatin, and pool it into a new 15-mL tube.
20. Prepare single-use 100-μL aliquots of the pooled supernatant containing the chromatin, and snap-freeze them in LN₂ for storage at -80°C.

Chromatin Quality and Quantity Determination and Input Chromatin Preparation

21. To determine the quality and quantity of the chromatin preparation, take a single 100-μL aliquot and reverse cross-link by incubation overnight at 65°C in a hybridization oven.
22. Clean up the test sample using DNA purification spin columns (QIAGEN catalog no. 28104).
Use only 50 μL per column because the amount of DNA may sometimes exceed the column capacity.
23. Collect the chromatin upon elution using 30 μL of Elution buffer, and repeat this step twice before combining both eluates into a single Eppendorf tube.
24. Mix 25 μL of the eluted chromatin with 5 μL of DNA loading buffer, and load 5, 10, and 15 μL of this mixture into the separate wells of a 1% agarose gel. Load sonicated salmon sperm DNA as a control to determine the amount of extracted DNA.
 - i. Mix 2.5 μL of a 1:10 diluted salmon sperm DNA stock of 1 mg/mL with 22.5 μL of ddH₂O and 5 μL of DNA loading dye.
 - ii. Load 5, 10, and 15 μL of the sonicated salmon sperm DNA mixture into separate wells of a 1% agarose gel.
25. Separate the sonicated chromatin and salmon sperm DNA by gel electrophoresis.
Chromatin DNA should appear as a smear ranging from 100 bp to 1 kb, with a maximum intensity around 500 bp.
26. If the chromatin has not been sufficiently sheared to the desired size, repeat the sample sonication as described in Step 17.
27. Quantify the reverse cross-linked and purified chromatin using a NanoDrop Spectrophotometer.

Blocking and Antibody Coating of Magnetic Beads

28. Add 1 mL of Block solution (PBS containing 0.5% BSA) to 100 μL of Dynabeads in a 1.5-mL Eppendorf tube.
29. Collect the Dynabeads using a magnetic rack; remove the supernatant, and wash the beads in 1.5 mL of Block solution for a total of three times.
30. Add 10 μg of antibody to the Dynabeads resuspended in 250 μL of Block solution and incubate overnight at 4°C.
31. Wash the Dynabeads three times in Block solution, and discard the Block solution.
32. Wash the blocked Dynabeads with RIPA buffer and aspirate. Store the Dynabeads in 100 μL of RIPA buffer at 4°C.

Setting Up the Chromatin Immunoprecipitation

33. Centrifuge the chromatin from Step 20 at 16,000g for 15 min at 4°C, and transfer it to a new tube.
34. Set up the chromatin immunoprecipitation as described in Protocol: **Immunoprecipitation** (DeCaprio and Kohl 2017), and incubate it overnight at 4°C.

Generally, 1 μg of antibody is added to 12.5 μg of chromatin DNA per immunoprecipitation, translating into the use of 10 μL of antibody prebound to Protein A magnetic beads.

Washing

35. Collect the Dynabeads using a magnetic rack, and wash the Dynabeads three times with 1 mL of RIPA buffer washes allowing for 1–2 min between washes.
36. Collect the Dynabeads using a magnetic rack, and wash them three times with 1 mL of EBC lysis buffer.
37. Collect the Dynabeads using a magnetic rack, and wash them twice with 1 mL of TE buffer containing 50 mM NaCl.
38. Collect the Dynabeads by centrifugation at 1000g for 3 min at 4°C, and pipette off any remaining TE buffer.

Dynabead collection using the magnetic rack does not ensure for complete removal of TE buffer stuck to the sides of the tube. Thus, centrifugation of the beads is recommended for complete buffer collection followed by removal.

Elution

39. Elute the sample from the Dynabeads by addition of 100 µL of T₅₀E₁₀S₁. Resuspend the beads by pipetting up and down about 10 times.

Be careful not to lose any beads in the pipette tip.

40. Incubate the resuspended Dynabeads for 30 min at 65°C while vortexing every 5–10 min.

Do not allow the beads to settle down.

41. Centrifuge the sample at 16,000g for 1 min at room temperature, and transfer the supernatant to a new 1.5-mL Eppendorf tube.

42. Elute the sample from the Dynabeads by addition of 100 µL of T₅₀E₁₀ (no SDS), and incubate the resuspended Dynabeads for 30 min at 65°C while vortexing every 5–10 min.

Do not allow the beads to settle down.

43. Centrifuge the sample at 16,000g for 1 min at room temperature, and transfer the supernatant to the Eppendorf tube containing the first elution (Step 41) to obtain a total volume of 200 µL.

Reverse Cross-Linking

44. Incubate the sample collected in Step 43 overnight at 65°C in a hybridization oven to reverse the cross-linking.
45. Add 20 µg of RNase A to the sample and incubate for 1 h at 37°C.
46. Add 40 µg of Proteinase K and incubate for 2 h at 65°C.

DNA Purification, Amplification, and Sequencing

47. Purify the sample using polymerase chain reaction (PCR) purification spin columns, and elute the DNA in 50 µL of T₅₀E₁₀ buffer. The eluted DNA can be analyzed by PCR or real-time PCR or stored at –80°C for subsequent use.

For Troubleshooting, see Table 1.

Dual Cross-Linking for Chromatin Immunoprecipitation

The dual cross-linking procedure serves the purpose of cross-linking proteins to proteins before the formaldehyde fixation of chromatin to protein and is particularly valuable for low-abundance protein antigens or when proteins interact indirectly with chromatin through associated proteins (Schmidt et al. 1977; Choe et al. 2002; Zeng et al. 2006; Adelmant and Marto 2009). This method, however, is of lesser value when protein targets are highly abundant, and, in this case, the main protocol is the better choice for a chromatin precipitation. The dual cross-linking procedure described here incorporates the use of disuccinimidyl glutarate [di(N-succinimidyl)] glutarate (DSG) during the initial cross-linking step before formaldehyde fixation.

TABLE 1. Troubleshooting

Problem	Possible cause	Solution
No specific signal observed	Inefficient precipitation of the cross-linked target protein	Increase the amount of antibody. Use a monoclonal antibody recognizing a different epitope. Use a different antibody with higher affinity to the target. Use an affinity-purified polyclonal antibody.
	Inefficient chromatin/protein cross-linking	Increase the cross-linking time. If protein/chromatin interaction is transient, explore the dual cross-linking method.
	Antibody epitope is masked by cross-linked chromatin.	Use a monoclonal antibody recognizing a different epitope. Use an affinity-purified polyclonal antibody.
	DNA fragments too large or too small	Monitor DNA fragment sizes on an agarose gel.
	Formaldehyde cross-linking time has been too short or too long.	Evaluate various cross-linking time periods.
	Wash buffer disrupted antibody–antigen interaction	Use milder Wash buffer. Reduce the salt and detergent concentration, or use a different detergent and reduce the number of washes.
High background in negative controls	Nonspecific sample binding to control antibodies	Clarify the lysate before addition of the specific antibody by centrifugation at 10,000g for 10 min at 4°C, transfer the supernatant to a clean tube without disturbing the pellet, and repeat the centrifugation.
	Genomic or plasmid contamination of reagents/buffers	Perform control PCRs without DNA templates.

For additional information, please refer to Table 1 in Protocol: **Immunoprecipitation** (DeCaprio and Kohl 2017).

48. Remove the growth medium from five 15-cm² cell culture dishes of proliferating cells by aspiration.

Avoid the use of overconfluent or dead cells.

49. Add 2 mL of PBS containing 1 mM PMSF to each cell culture dish; scrape and transfer the cells to a 50-mL conical tube before cell pelleting by centrifugation at 400g for 5 min at room temperature.
50. Wash the pellet in 20 mL of PBS with 1 mM PMSF and centrifuge at 400g for 5 min at room temperature. Repeat this step for a total of three washes. Resuspend the cell pellet in 20 mL of PBS with 1 mM PMSF by gentle pipetting.

Ensure that no cell clumps are present and that all cells have been properly resuspended.

51. Add 80 µL of 0.5 M disuccinimidyl glutarate (DSG) to the 20 mL of PBS/PMSF-resuspended cells, resulting in a final DSG concentration of 2 mM, and mix immediately.
52. Incubate the sample for 45 min at room temperature on a rotating wheel set to medium speed (5–10 rpm).
53. Collect the cells by centrifugation at 400g for 10 min at room temperature.
54. Wash the cell pellet in 20 mL of PBS, and centrifuge it at 400g for 5 min at room temperature. Repeat this step for a total of three washes.

It is extremely important to wash off and remove the DSG completely.

55. Resuspend the cell pellet in 20 mL of PBS at room temperature.
56. Add 540 µL of 37% formaldehyde. Mix immediately and incubate the sample for 15 min at room temperature on a rotating wheel set to medium speed (5–10 rpm).
57. Quench the reaction by addition of 1 mL of 2.5 M glycine. Mix immediately, and incubate the sample for 10 min at room temperature on a rotating wheel set to medium speed (5–10 rpm).
58. Collect the cells by centrifugation at 400g for 5 min at 4°C, and aspirate the supernatant.
59. Wash the cell pellet twice with 10 mL of ice-cold PBS while keeping the samples on ice.
60. Snap-freeze the cells in LN₂ for storage at –80°C.

61. Proceed with the chromatin preparation and immunoprecipitation as described above, starting at Step 14.

For Troubleshooting, see Table 1.

DISCUSSION

The chromatin immunoprecipitation technique represents the most commonly applied method for investigating specific interactions between proteins and DNA on a genome-wide scale. Developed and refined in the mid-1980s, this technique is capable of elucidating the various underlying mechanisms of gene regulation, replication, and DNA damage repair (Takahashi et al. 2000; Wang et al. 2002; Weinmann and Farnham 2002; Aparicio et al. 2004; Viens et al. 2004; Breiling and Orlando 2006; Havis et al. 2006; Jin et al. 2006; Nelson et al. 2006; O'Neill et al. 2006; Sandmann et al. 2006; Haring et al. 2007; Tansey 2007; Brand et al. 2008; Cho et al. 2008; Dahl and Collas 2008a,b; Gossett and Lieb 2008; Jiang et al. 2008; Matevossian and Akbarian 2008; Mukhopadhyay et al. 2008; Saleh et al. 2008; Bermejo et al. 2009; Caby and Pierce 2009; Fode and Gatz 2009; Furlan-Magaril et al. 2009; Gaillard et al. 2009; Lindeman et al. 2009; Magdalena and Goval 2009; Papantonis and Lecanidou 2009; Percipalle and Obrdlik 2009; Sikes et al. 2009; Sorensen and Collas 2009; Tong and Falk 2009; Buro et al. 2010; Grably and Engelberg 2010; Haudry et al. 2010; Jayani et al. 2010; Kaufmann et al. 2010; Pien and Grossniklaus 2010; Ricardi et al. 2010; Villar and Kohler 2010; Adli and Bernstein 2011; Chernukhin et al. 2011; Collas 2011; Cosseau and Grunau 2011; de Folter 2011; de Medeiros 2011; Fanelli et al. 2011; Floriddia et al. 2011; Kulyassov et al. 2011; Strenkert et al. 2011; Akkers et al. 2012; Bonn et al. 2012; Cheung et al. 2012; Duc et al. 2012; Gilfillan et al. 2012; Heliot and Cereghini 2012; Jung and Adjaye 2012; Lin et al. 2012; Micsinai et al. 2012; Murgatroyd et al. 2012; Ramaswami et al. 2012; Sailaja et al. 2012; Schuch et al. 2012; Singh and Szabo 2012; Wu et al. 2012; Adhikary and Muller 2013; Bolli et al. 2013; Blecher-Gonen et al. 2013; Cheaib and Simon 2013; Le May and Li 2013; Liao and Lu 2013; Lopez-Rubio et al. 2013; Peng and Chen 2013; van den Boogaard et al. 2013).

Numerous ChIP-validated polyclonal and monoclonal antibodies have been developed specific to the epigenetic (“control from above”) regulation of gene expression and DNA replication that target DNA and histone modifications, including methylation, acetylation, ubiquitylation, phosphorylation, and sumoylation. Although the quality and specificity of these and other antibodies as well as other aspects of the immunoprecipitation process are important to the success of a chromatin immunoprecipitation, the quality of the sheared chromatin material extracted from cells or tissues is essential to the overall goal of a chromatin immunoprecipitation. Various methods for the effective shearing of chromatin have been developed and include the traditional sonication and enzymatic digestion approaches and the recently developed controllable isothermal shearing. The latter overcomes the shortcomings of sonication or enzymatic digestion by foregoing epitope damage using excessive heating and the introduction of shearing biases, thus consistently delivering reproducible results.

Although ChIP provides valuable information, the limitations of the assay need to be taken into consideration. Polyclonal or monoclonal antibodies used for the immunoprecipitation of specific target proteins might cross-react with other nuclear proteins, resulting in the misinterpretation of the results. The ChIP assay is not a functional assay and does not elucidate the function of the protein associated with the immunoprecipitated DNA, nor can it reveal the exact binding site of the protein during its interaction with the DNA. Furthermore, ChIP assays often deliver low signals in comparison with the negative controls and could result in possible misinterpretation.

Nonetheless, since its development, refinement, and modification, ChIP technology has resulted in the development of two further assays, commonly referred to as ChIP-chip (or ChIP-on-chip) and ChIP-Seq. The ChIP-chip assay requires the amplification of purified DNA fragments followed by the annealing of fragments to DNA chips or microarrays composed of immobilized DNA oligonucleotides covering the entire genome, thus providing valuable information about the genomic region(s) specifically interacting with the immunoprecipitated protein. Unlike the ChIP-chip assay, the ChIP-sequence (ChIP-Seq) assay is not influenced by variable or weak DNA annealing efficiencies and is

considered to be more accurate by providing quantitative data, including the actual DNA sequences specifically interacting with the immunoprecipitated protein. Together, all three assays—ChIP, ChIP-chip, and ChIP-Seq—can provide valuable data and information about the proteins associated with DNA and their function in gene regulation.

RECIPES

Disuccinimidyl Glutarate (DSG; 0.5 M)

Reagent	Amount per 300 μ L of solution (w/v)	Final concentration
Disuccinimidyl glutarate	50 mg	0.5 M
DMSO	300 μ L	

Store at -20°C .

DNA Loading Buffer (6 \times)

30% (v/v) glycerol
0.25% (w/v) bromophenol blue
0.25% (w/v) xylene cyanol FF

Store at 4°C .

EBC Lysis Buffer for ChIP

Reagent	Volume per 100 mL of solution (v/v)	Final concentration
NaCl (5 M)	2.4 mL	120 mM
Nonidet P-40 (10%)	5.0 mL	0.5%
Leupeptin (10 mg/mL)	50 μ L	5 μ g/mL
Aprotinin (10 mg/mL)	100 μ L	10 μ g/mL
PMSF (0.1 M)	1 mL	1 mM
Na-vanadate (0.2 M)	100 μ L	0.2 mM
NaF (1 M)	10 mL	100 mM
Tris-HCl (1 M, pH 8.0)	5.0 mL	50 mM
ddH ₂ O	Add to 100 mL	

Store at 4°C . Add protease inhibitors immediately before use.

Elution Buffer for ChIP

50 mM Tris-Cl (pH 8.0)
10 mM EDTA
1% SDS

Store at room temperature.

Formaldehyde Solution (11%)

Reagent	Volume per 50 mL of solution (v/v)	Final concentration
HEPES-KOH (1 M, pH 7.5)	2.5 mL	50 mM
NaCl (5 M)	1.0 mL	100 mM
EDTA (0.5 M)	100 μ L	1 mM
EGTA (0.5 M)	50 μ L	0.5 mM
Formaldehyde (37%)	14.9 mL	11%
ddH ₂ O	Add to 100 mL	

Store at room temperature.

Lysis Buffer 1 for ChIP

Reagent	Volume per 100 mL of solution (v/v)	Final concentration
HEPES-KOH (1 M, pH 7.5)	5 mL	50 mM
NaCl (5 M)	2.8 mL	140 mM
EDTA (0.5 M)	0.2 mL	1 mM
Glycerol (50%)	20 mL	10%
Nonidet P-40 (10%)	5.0 mL	0.5%
Triton X-100 (10%)	2.5 mL	0.25%
ddH ₂ O	Add to 100 mL	
Store at 4°C.		

Lysis Buffer 2 for ChIP

Reagent	Volume per 100 mL of solution (v/v)	Final concentration
Tris-HCl (1 M, pH 8.0)	1.0 mL	10 mM
NaCl (5 M)	4.0 mL	200 mM
EDTA (0.5 M)	0.2 mL	1 mM
EGTA (0.5 M)	0.1 mL	0.5 mM
ddH ₂ O	Add to 100 mL	
Store at 4°C.		

Lysis Buffer 3 for ChIP

Reagent	Volume per 100 mL of solution (v/v)	Final concentration
Tris-HCl (1 M, pH 8.0)	1.0 mL	10 mM
NaCl (5 M)	2.0 mL	100 mM
EDTA (0.5 M)	0.2 mL	1 mM
EGTA (0.5 M)	0.1 mL	0.5 mM
Na-deoxycholate (10%)	1.0 mL	0.1%
N-Lauroyl-sarcosine (20%)	2.5 mL	0.5%
ddH ₂ O	Add to 100 mL	
Store at 4°C.		

RIPA Buffer Solution for ChIP

Reagent	Volume per 500 mL of solution (v/v)	Final concentration
HEPES-KOH (1 M, pH 7.6)	25 mL	50 mM
LiCl (5 M)	50 mL	500 mM
EDTA (0.5 M)	1 mL	1 mM
Nonidet P-40 (10%)	50 mL	1%
Na-deoxycholate (10%)	35 mL	0.7%
ddH ₂ O	Add to 500 mL	

Prepared RIPA buffer should be aliquoted and stored at −20°C. Add protease and/or phosphatase inhibitors to a thawed aliquot before immediate use. Discard and do not freeze again.

REFERENCES

- Adelmant G, Marto JA. 2009. Protein complexes: The forest and the trees. *Expert Rev Proteomics* 6: 5–10.
- Adhikary T, Muller R. 2013. In vivo studies of PPAR–chromatin interactions: Chromatin immunoprecipitation for single-locus and genome-wide analyses. *Methods Mol Biol* 952: 175–185.
- Adli M, Bernstein BE. 2011. Whole-genome chromatin profiling from limited numbers of cells using nano-ChIP-seq. *Nat Protoc* 6: 1656–1668.
- Akkers RC, Jacobi UG, Veenstra GJ. 2012. Chromatin immunoprecipitation analysis of *Xenopus* embryos. *Methods Mol Biol* 917: 279–292.

- Aparicio O, Geisberg JV, Struhl K. 2004. Chromatin immunoprecipitation for determining the association of proteins with specific genomic sequences in vivo. *Curr Protoc Cell Biol* 23: 17.7.1–17.7.23.
- Bermejo R, Katou YM, Shirahige K, Foiani M. 2009. ChIP-on-chip analysis of DNA topoisomerases. *Methods Mol Biol* 582: 103–118.
- Blecher-Gonen R, Barnett-Itzhaki Z, Jaitin D, Amann-Zalcenstein D, Lara-Astiaso D, Amit I. 2013. High-throughput chromatin immunoprecipitation for genome-wide mapping of in vivo protein–DNA interactions and epigenomic states. *Nat Protoc* 8: 539–554.
- Bolli P, Vardabasso C, Bernstein E, Chaudhry HW. 2013. Chromatin immunoprecipitation of adult murine cardiomyocytes. *Curr Protoc Cell Biol* 58: 17.14.1–17.14.16.
- Bonn S, Zinzen RP, Perez-Gonzalez A, Riddell A, Gavin AC, Furlong EE. 2012. Cell type-specific chromatin immunoprecipitation from multicellular complex samples using BiTS-ChIP. *Nat Protoc* 7: 978–994.
- Brand M, Rampalli S, Chaturvedi CP, Dilworth FJ. 2008. Analysis of epigenetic modifications of chromatin at specific gene loci by native chromatin immunoprecipitation of nucleosomes isolated using hydroxyapatite chromatography. *Nat Protoc* 3: 398–409.
- Breiling A, Orlando V. 2006. Chromatin immunoprecipitation (ChIP) of protein complexes: Mapping of genomic targets of nuclear proteins in cultured cells. *Cold Spring Harb Protoc* doi: 10.1101/pdb.prot4560.
- Buro LJ, Shah S, Henriksen MA. 2010. Chromatin immunoprecipitation (ChIP) to assay dynamic histone modification in activated gene expression in human cells. *J Vis Exp* 2010: 2053.
- Caby S, Pierce RJ. 2009. Quantitative chromatin immunoprecipitation (Q-ChIP) applied to *Schistosoma mansoni*. *Mol Biochem Parasitol* 166: 77–80.
- Cheab M, Simon M. 2013. Dynamic chromatin remodelling of ciliate macronuclear DNA as determined by an optimized chromatin immunoprecipitation (ChIP) method for *Paramecium tetraurelia*. *Appl Microbiol Biotechnol* 97: 2661–2670.
- Chernukhin I, Kang SY, Brown S, Gretton S, Mendez-Catala CF, Cowieson D, Klenova E. 2011. BioVyon Protein A, an alternative solid-phase affinity matrix for chromatin immunoprecipitation. *Anal Biochem* 412: 183–188.
- Cheung HH, Lee TL, Rennert OM, Chan WY. 2012. Methylation profiling using methylated DNA immunoprecipitation and tiling array hybridization. *Methods Mol Biol* 825: 115–126.
- Cho BK, Knight EM, Palsson BO. 2008. Genomewide identification of protein binding locations using chromatin immunoprecipitation coupled with microarray. *Methods Mol Biol* 439: 131–145.
- Choe W, Budd M, Imamura O, Hoopes L, Campbell JL. 2002. Dynamic localization of an Okazaki fragment processing protein suggests a novel role in telomere replication. *Mol Cell Biol* 22: 4202–4217.
- Collas P. 2011. A chromatin immunoprecipitation protocol for small cell numbers. *Methods Mol Biol* 791: 179–193.
- Cosseau C, Grunau C. 2011. Native chromatin immunoprecipitation. *Methods Mol Biol* 791: 195–212.
- Dahl JA, Collas P. 2008a. MicroChIP—A rapid micro chromatin immunoprecipitation assay for small cell samples and biopsies. *Nucleic Acids Res* 36: 17.
- Dahl JA, Collas P. 2008b. A rapid micro chromatin immunoprecipitation assay (microChIP). *Nat Protoc* 3: 1032–1045.
- DeCaprio J, Kohl TO. 2017. Immunoprecipitation. *Cold Spring Harb Protoc* doi: 10.1101/pdb.prot098640.
- de Folter S. 2011. Protein tagging for chromatin immunoprecipitation from *Arabidopsis*. *Methods Mol Biol* 678: 199–210.
- de Medeiros RB. 2011. Sequential chromatin immunoprecipitation assay and analysis. *Methods Mol Biol* 791: 225–237.
- Duc TN, Hassanzadeh-Ghassabeh G, Saerens D, Peeters E, Charlier D, Muyldermans S. 2012. Nanobody-based chromatin immunoprecipitation. *Methods Mol Biol* 911: 491–505.
- Fanelli M, Amatori S, Barozzi I, Minucci S. 2011. Chromatin immunoprecipitation and high-throughput sequencing from paraffin-embedded pathology tissue. *Nat Protoc* 6: 1905–1919.
- Floriddia E, Nguyen T, Di Giovanni S. 2011. Chromatin immunoprecipitation from dorsal root ganglia tissue following axonal injury. *J Vis Exp* 2011: 2803.
- Fode B, Gatz C. 2009. Chromatin immunoprecipitation experiments to investigate in vivo binding of *Arabidopsis* transcription factors to target sequences. *Methods Mol Biol* 479: 261–272.
- Furlan-Magaril M, Rincon-Arango H, Recillas-Targa F. 2009. Sequential chromatin immunoprecipitation protocol: ChIP-reChIP. *Methods Mol Biol* 543: 253–266.
- Gaillard H, Wellinger RE, Aguilera A. 2009. Methods to study transcription-coupled repair in chromatin. *Methods Mol Biol* 523: 141–159.
- Gilfillan GD, Hughes T, Sheng Y, Hjorthaug HS, Straub T, Gervin K, Harris JR, Undlien DE, Lyle R. 2012. Limitations and possibilities of low cell number ChIP-seq. *BMC Genomics* 13: 1471–2164.
- Gossett AJ, Lieb JD. 2008. DNA Immunoprecipitation (DIP) for the determination of DNA-binding specificity. *Cold Spring Harb Protoc* doi: 10.1101/pdb.prot4972.
- Grably M, Engelberg D. 2010. A detailed protocol for chromatin immunoprecipitation in the yeast *Saccharomyces cerevisiae*. *Methods Mol Biol* 638: 211–224.
- Haring M, Offermann S, Danker T, Horst I, Peterhansel C, Stam M. 2007. Chromatin immunoprecipitation: Optimization, quantitative analysis and data normalization. *Plant Methods* 3: 11.
- Haudry Y, Ramialison M, Paten B, Wittbrodt J, Ettwiller L. 2010. Using Trawler_standalone to discover overrepresented motifs in DNA and RNA sequences derived from various experiments including chromatin immunoprecipitation. *Nat Protoc* 5: 323–334.
- Havis E, Anselme I, Schneider-Maunoury S. 2006. Whole embryo chromatin immunoprecipitation protocol for the in vivo study of zebrafish development. *BioTechniques* 40: 34–40.
- Heliot C, Cereghini S. 2012. Analysis of in vivo transcription factor recruitment by chromatin immunoprecipitation of mouse embryonic kidney. *Methods Mol Biol* 886: 275–291.
- Jayani RS, Ramanujam PL, Galande S. 2010. Studying histone modifications and their genomic functions by employing chromatin immunoprecipitation and immunoblotting. *Methods Cell Biol* 98: 35–56.
- Jiang Y, Matevosian A, Huang HS, Straubhaar J, Akbarian S. 2008. Isolation of neuronal chromatin from brain tissue. *BMC Neurosci* 9: 1471–2202.
- Jin VX, Rabinovich A, Squazzo SL, Green R, Farnham PJ. 2006. A computational genomics approach to identify cis-regulatory modules from chromatin immunoprecipitation microarray data—A case study using E2F1. *Genome Res* 16: 1585–1595.
- Jung M, Adjaye J. 2012. Chromatin immunoprecipitation-based analysis of gene regulatory networks operative in human embryonic stem cells. *Methods Mol Biol* 873: 269–280.
- Kaufmann K, Muino JM, Osteras M, Farinelli L, Krajewski P, Angenent GC. 2010. Chromatin immunoprecipitation (ChIP) of plant transcription factors followed by sequencing (ChIP-SEQ) or hybridization to whole genome arrays (ChIP-CHIP). *Nat Protoc* 5: 457–472.
- Kulyassov A, Shoaib M, Ogryzko V. 2011. Use of in vivo biotinylation for chromatin immunoprecipitation. *Curr Protoc Cell Biol* 51: 17.12.1–17.12.22.
- Le May M, Li Q. 2013. Analysis of p300 occupancy at the early stage of stem cell differentiation by chromatin immunoprecipitation. *Methods Mol Biol* 977: 315–321.
- Liao JM, Lu H. 2013. ChIP for identification of p53 responsive DNA promoters. *Methods Mol Biol* 962: 201–210.
- Lin X, Tirichine L, Bowler C. 2012. Protocol: Chromatin immunoprecipitation (ChIP) methodology to investigate histone modifications in two model diatom species. *Plant Methods* 8: 48.
- Lindeman LC, Vogt-Kielland LT, Alestrom P, Collas P. 2009. Fish'n ChIPs: Chromatin immunoprecipitation in the zebrafish embryo. *Methods Mol Biol* 567: 75–86.
- Lopez-Rubio JJ, Siegel TN, Scherf A. 2013. Genome-wide chromatin immunoprecipitation-sequencing in *Plasmodium*. *Methods Mol Biol* 923: 321–333.
- Magdalena J, Goyal JJ. 2009. Methyl DNA immunoprecipitation. *Methods Mol Biol* 567: 237–247.
- Matevosian A, Akbarian S. 2008. A chromatin assay for human brain tissue. *J Vis Exp* 21: 717.
- Micsinai M, Parisi F, Strino F, Asp P, Dynlacht BD, Kluger Y. 2012. Picking ChIP-seq peak detectors for analyzing chromatin modification experiments. *Nucleic Acids Res* 40: 3.
- Mukhopadhyay A, Deplancke B, Walhout AJ, Tissenbaum HA. 2008. Chromatin immunoprecipitation (ChIP) coupled to detection by quantitative real-time PCR to study transcription factor binding to DNA in *Caenorhabditis elegans*. *Nat Protoc* 3: 698–709.

J. DeCaprio and T.O. Kohl

- Murgatroyd C, Hoffmann A, Spengler D. 2012. In vivo ChIP for the analysis of microdissected tissue samples. *Methods Mol Biol* **809**: 135–148.
- Nelson JD, Denisenko O, Bomsztyk K. 2006. Protocol for the fast chromatin immunoprecipitation (ChIP) method. *Nat Protoc* **1**: 179–185.
- O'Neill LP, VerMilyea MD, Turner BM. 2006. Epigenetic characterization of the early embryo with a chromatin immunoprecipitation protocol applicable to small cell populations. *Nat Genet* **38**: 835–841.
- Papantonis A, Lecanidou R. 2009. A modified chromatin-immunoprecipitation protocol for silkworm ovarian follicular cells reveals C/EBP and GATA binding modes on an early chorion gene promoter. *Mol Biol Rep* **36**: 733–736.
- Peng GH, Chen S. 2013. Double chromatin immunoprecipitation: Analysis of target co-occupancy of retinal transcription factors. *Methods Mol Biol* **935**: 311–328.
- Percipalle P, Obrdlik A. 2009. Analysis of nascent RNA transcripts by chromatin RNA immunoprecipitation. *Methods Mol Biol* **567**: 215–235.
- Pien S, Grossniklaus U. 2010. Chromatin immunoprecipitation protocol for histone modifications and protein–DNA binding analyses in *Arabidopsis*. *Methods Mol Biol* **631**: 209–220.
- Ramaswami S, Manna S, Juvekar A, Kennedy S, Vancura A, Vancurova I. 2012. Chromatin immunoprecipitation analysis of NFκB transcriptional regulation by nuclear IκBα in human macrophages. *Methods Mol Biol* **809**: 121–134.
- Ricardi MM, Gonzalez RM, Iusem ND. 2010. Protocol: Fine-tuning of a chromatin immunoprecipitation (ChIP) protocol in tomato. *Plant Methods* **6**: 11.
- Sailaja BS, Takizawa T, Meshorer E. 2012. Chromatin immunoprecipitation in mouse hippocampal cells and tissues. *Methods Mol Biol* **809**: 353–364.
- Saleh A, Alvarez-Venegas R, Avramova Z. 2008. An efficient chromatin immunoprecipitation (ChIP) protocol for studying histone modifications in *Arabidopsis* plants. *Nat Protoc* **3**: 1018–1025.
- Sandmann T, Jakobsen JS, Furlong EE. 2006. ChIP-on-chip protocol for genome-wide analysis of transcription factor binding in *Drosophila melanogaster* embryos. *Nat Protoc* **1**: 2839–2855.
- Schmidt HJ, Schaum U, Pichotka JP. 1977. [The influence of mode and intensity of homogenization on the absolute value and stability of oxygen consumption of guinea pig liver homogenates (author's transl)]. *Z Naturforsch C* **32**: 908–912.
- Schuch R, Agelopoulos K, Neumann A, Brandt B, Burger H, Korsching E. 2012. Site-specific chromatin immunoprecipitation: A selective method to individually analyze neighboring transcription factor binding sites in vivo. *BMC Res Notes* **5**: 109.
- Sikes ML, Bradshaw JM, Ivory WT, Lunsford JL, McMillan RE, Morrison CR. 2009. A streamlined method for rapid and sensitive chromatin immunoprecipitation. *J Immunol Methods* **344**: 58–63.
- Singh P, Szabo PE. 2012. Chromatin immunoprecipitation to characterize the epigenetic profiles of imprinted domains. *Methods Mol Biol* **925**: 159–172.
- Sorensen AL, Collas P. 2009. Immunoprecipitation of methylated DNA. *Methods Mol Biol* **567**: 249–262.
- Strenkert D, Schmollinger S, Schroda M. 2011. Protocol: Methodology for chromatin immunoprecipitation (ChIP) in *Chlamydomonas reinhardtii*. *Plant Methods* **7**: 35.
- Takahashi K, Saitoh S, Yanagida M. 2000. Application of the chromatin immunoprecipitation method to identify in vivo protein–DNA associations in fission yeast. *Sci STKE* **2000**: p11.
- Tansey WP. 2007. Yeast chromatin immunoprecipitation (ChIP) assay. *Cold Spring Harb Protoc* doi: 10.1101/pdb.prot4642.
- Tong Y, Falk J. 2009. Genome-wide analysis for protein–DNA interaction: ChIP-chip. *Methods Mol Biol* **590**: 235–251.
- van den Boogaard M, Wong LY, Christoffels VM, Barnett P. 2013. Acquisition of high quality DNA for massive parallel sequencing by in vivo chromatin immunoprecipitation. *Methods Mol Biol* **977**: 53–64.
- Viens A, Mechold U, Lehrmann H, Harel-Bellan A, Ogryzko V. 2004. Use of protein biotinylation in vivo for chromatin immunoprecipitation. *Anal Biochem* **325**: 68–76.
- Villar CB, Kohler C. 2010. Plant chromatin immunoprecipitation. *Methods Mol Biol* **655**: 401–411.
- Wang H, Tang W, Zhu C, Perry SE. 2002. A chromatin immunoprecipitation (ChIP) approach to isolate genes regulated by AGL15, a MADS domain protein that preferentially accumulates in embryos. *Plant J* **32**: 831–843.
- Weinmann AS, Farnham PJ. 2002. Identification of unknown target genes of human transcription factors using chromatin immunoprecipitation. *Methods* **26**: 37–47.
- Wu AR, Kawahara TL, Rapicavoli NA, van Riggelen J, Shroff EH, Xu L, Felsher DW, Chang HY, Quake SR. 2012. High throughput automated chromatin immunoprecipitation as a platform for drug screening and antibody validation. *Lab Chip* **12**: 2190–2198.
- Zeng PY, Vakoc CR, Chen ZC, Blobel GA, Berger SL. 2006. In vivo dual cross-linking for identification of indirect DNA-associated proteins by chromatin immunoprecipitation. *Biotechniques* **41**: 694–698.



Cold Spring Harbor Protocols

Chromatin Immunoprecipitation

James DeCaprio and Thomas O. Kohl

Cold Spring Harb Protoc; doi: 10.1101/pdb.prot098665

Email Alerting Service

Receive free email alerts when new articles cite this article - [click here](#).

Subject Categories

Browse articles on similar topics from *Cold Spring Harbor Protocols*.

[DNA:Protein Interactions](#) (74 articles)

[Immunoprecipitation](#) (75 articles)

[Molecular Biology, general](#) (1314 articles)

To subscribe to *Cold Spring Harbor Protocols* go to:
<http://cshprotocols.cshlp.org/subscriptions>
