

A Quick, No Frills Approach to Mouse Genotyping

Manuel E. Lopez*

Department of Developmental Biology, Stanford University, Stanford, USA

*For correspondence: manduardo@gmail.com

[Abstract] Mice are extremely powerful mammalian genetic model organisms for basic and medical research, but managing a colony of transgenic mice is time consuming and expensive, many times requiring the help of dedicated technicians. Slow and laborious genotyping procedures add to the hassle. Outsourcing is costly and may not be as fast as desired, especially when setting up time sensitive experiments. Ultrafast genotyping protocols often require real-time PCR instruments and commercial reagents that may not be economical or practical. This protocol, adapted from methods suggested by The Jackson Laboratory, employs a minimalist approach that maximizes convenience by simplifying the tissue digestion/DNA extraction process and using a high-speed electrophoresis system for sample analysis. Genotype PCR results can be obtained in 3 h or less for as many samples as can fit in a PCR machine or can be efficiently handled by a user. Subsequent ethanol or chloroform purified DNA can be used in a standard PCR reaction to roughly identify a homozygous and a hemizygous mouse.

Materials and Reagents

1. NaOH (NaOH pellets)
2. Taq DNA Polymerase with ThermoPol buffer (New England Biolabs, catalog number: M0267X*, M0267L, or M0267S)

**Note: At 4,000 cycles, Taq DNA Polymerase becomes a limiting reagent. ThermoPol buffer recipe is available at NEB website. This buffer can be ordered separately from NEB (New England Biolabs, catalog number: B9004S)*

3. Primers, recommend to be 18-21 bp in length, have a melting temperature above 56 °C and around 58 °C, and produce amplicons of 150-600 bp.
4. DNA loading buffer
Recommend Orange G (Sigma-Aldrich, catalog number: O3756) instead of bromophenol blue for loading dye
5. DNA ladder range for 100-800 bp range
Recommend 1 kb Plus DNA ladder (Life Technologies, catalog number: 10787-018)
6. Agarose

7. EtBr (Sigma-Aldrich, catalog number: E8751) or SYBR Safe DNA Gel Stain (Life Technologies, catalog number: S33102)
8. NaOAc
9. EtOH
10. Phenol/Chloroform/Isoamyl alcohol (25: 24: 1) (Do not use acid Phenol)
11. 50 mM NaOH in dH₂O (see Recipes)
12. 10 mM dNTP Mix (see Recipes)
13. 1x TAE buffer (see Recipes)
14. dNTP set 100 mM each A, C, G, T (GE Life Sciences, catalog number: 28406552) (see Recipes)
15. 0.3 M NaOAc in ddH₂O (see Recipes)

Equipment

1. PCR Thermal cycler (96 well capacity preferred)
2. Centrifuges (mini for PCR tubes and microcentrifuge for 1.5 ml tubes)
3. Liberty 1 buffer-less high speed gel system (Neuvitro, 6Mgel - SYS-LBT1)
4. Liberty 1, 12-channel pipette compatible 13 teeth combs (Neuvitro, 6Mgel - CMS-1315)
5. Multichannel Pipette 2-GE |
6. Repeat Pipettor 10-FG |
7. PCR tube with cap, 8 or 12 PCR strip tubes, or 96-well PCR plate

Procedure

A. Part I. Digest tissue for genotyping

Note: This protocol is performed on 1-48 samples at the same time using 12-strip PCR tubes. The strip PCR tubes allow the use of multichannel pipettes to transfer solutions. 96 PCR plates can be used as well.

1. Tissue digestion:
 - a. Place a roughly 2-3 mm ear clip, tail, or other tissue biopsy in each PCR tube as they are obtained.
 - b. *Caution: 50 mM NaOH is a strong base and can be corrosive. It is more difficult to use a repeat pipettor to squirt out a sufficient volume of 50 mM NaOH into each PCR tube.*

- c. Incubate in thermo cycler for 95 °C for 30 min to 1 h. 45 min is recommended.
Note: Samples have been heated for up to 2 h and as little as 15 min. 15-30 min is sufficient for neonatal tissue.
- d. Immediately flick tubes.
Key feature: Tissue sample should partially break apart when PCR tube is flicked several times while still hot. Liquid should be cloudy. If tissue is not falling apart, incubate for longer.
- e. Allow samples to cool to room temperature (RT).
- f. Briefly spin to remove liquid from caps.
- g. Proceed to Part II or store samples at room temperature overnight for next day use.
Notes:
 - i. *Samples can be kept at room temperature for up to 3 days. Samples left for over a week have been used successfully but not recommended.*
 - ii. *Freeze at -20 °C for long-term storage. Neutralization step for base pH created by NaOH is not required but may benefit long-term storage. 10x ThermoPol buffer or 1 volume of 0.3 M NaOAc in dH₂O can be used to neutralize.*

B. Part II. PCR reaction and gel electrophoresis

2. PCR after tissue digestion:

h. PCR master mix setup.

Note: Always make for more reactions than needed. For every 12 samples make for 2 samples extra, e.g. for 48 samples, make master mix for 56 reactions total.

Per reaction:

10 µl dH₂O

10 µl 10x ThermoPol buffer

10 µl 10 mM dNTPs

1 µl 10 mM MgCl₂

1 µl Taq DNA Polymerase

1 µl

- i. 95 °C for 30 sec
- j. 95 °C for 30 sec
- k. Briefly spin.
- l. Samples can be kept at 4 °C until a PCR machine is available same day.
- m. PCR program
 - i. 95 - 2 min
 - ii. 95 - 30 sec

Acknowledgments

This protocol was adapted from previous work (Lopez *et al.*, 2011; Lopez *et al.*, 2012).

References

1. Lopez, M. E., Klein, A. D., Dimbil, U. J. and Scott, M. P. (2011). [Anatomically defined neuron-based rescue of neurodegenerative Niemann-Pick type C disorder.](#) *J Neurosci* 31(12): 4367-4378.
2. Lopez, M. E., Klein, A. D., Hong, J., Dimbil, U. J. and Scott, M. P. (2012). [Neuronal and epithelial cell rescue resolves chronic systemic inflammation in the lipid storage disorder Niemann-Pick C.](#) *Hum Mol Genet* 21(13): 2946-2960.